



## Scaling tip for scaler, KaVo Piezoled 201 compatible

Scaling tip for supragingival tartar removal in all quadrants.  
Compatible KaVo Piezoled 201

**Available**

**SKU:** PG021K

**Categories:** KaVo

### Description

---

The scaling tip for scaler compatible with KaVo Piezoled 201 is designed for supragingival calculus removal in all quadrants. It is intended for professional scaling procedures, supporting the removal of calculus deposits and stains from supragingival tooth surfaces.

### Specifications

- Compatible with KaVo Piezoled 201
- For professional dental use