



Scaling tip for scaler, KaVo Piezoled 202 compatible

Scaling tip for supragingival and subgingival tartar removal in all quadrants. Compatible KaVo Piezoled 202

Available

SKU: PG022K

Categories: KaVo

Description

The scaling tip for scaler compatible with KaVo Piezoled 202 is a tip for dental scaling, indicated for removing calculus deposits and supragingival, subgingival and interdental stains. The tip is suitable for the daily work of dentists and hygienists, allowing access even to anatomically difficult areas while maintaining effective reach of the treatment area. Indicated for the removal of supragingival and subgingival calculus in all quadrants.

Specifications

- Compatible with KaVo Piezoled 202
- Suitable for use with KaVo Piezoled piezoelectric scalers
- For professional dental use