



## Sealing machine Termoseal 300, sealing 36 cm, continuous heating

Thermoseal 300 sealing machine, sealing length 36 cm, width 1.5 cm, continuous heating

## **Available**

SKU: SEA300

Categories: Sealing Machines

## **Description**

Termoseal 300 heat sealer - Continuous sealing and constant temperature for sterile pouches.

The **Termoseal 300 heat sealer** is a compact, reliable device for sealing sterile pouches in medical and dental settings. Its **continuous heating system** ensures consistent, uninterrupted sealing for maximum safety and workflow efficiency.

Precision and consistent control.

With a **constant temperature system**, Termoseal 300 delivers perfect results on all heat-sealable materials, preventing overheating and ensuring airtight, durable seals. The **exposed aluminium bar** is robust and virtually indestructible, ensuring long-term stability.

Functional and practical design.

Featuring a **sealing width of 15mm** and a **sealing length of 360mm**, Termoseal 300 includes a **375mm roll holder rod**for quick and easy material handling.

Instant readiness and continuous operation.

Its continuous cycle eliminates downtime between seals, making it ideal for clinics and laboratories with high sterilisation throughput.

## Supplied heat sealer.

• Thermal sealing machine



- n.2 Support rodsn.1 Roll support rod
- Power cable

| Features                | Code SEA300      |
|-------------------------|------------------|
| Weld width              | 15 mm            |
| Weld length             | 360 mm           |
| Roll support rod length | 375 mm           |
| Power                   | 450 watt         |
| Power supply            | 230 vac          |
| Measurements (LxDxH)    | 500X260X200/260h |
| Weight                  | 6.5 Kg           |
| Warranty (months)       | 24               |
| Certifications          | CE               |