



## Service unit 500X475X810h, 1 drawer and 1 door, Ral color handles

**Available**

**SKU:** DE9.32.1SO

**Categories:** Side Cabinets

### Description

---

The **modular clinical side cabinet with drawer and door**, measuring 500x475x810h, is designed to improve organisation within clinical spaces. Its compact form helps make efficient use of the available room, providing order, practicality and reliable support during daily activities.

### STRUCTURE AND FUNCTIONALITY

The cabinet features a clean and simple design, ideal for settings where hygiene and clarity are essential. The drawer and door allow tools and materials to be arranged logically, ensuring they are always easy to reach. The internal layout streamlines workflows and reduces preparation times during clinical procedures.

### CLINICAL MODULARITY

This clinical side cabinet with drawer and door fits seamlessly into different medical and multidisciplinary practices. Its modular nature allows it to adapt to various clinical layouts, improving efficiency and helping maintain a professional and welcoming environment for both clinicians and patients.

### MOBILITY AND OPERATIONAL SUPPORT

Fitted with wheels, the cabinet can be moved quickly and positioned where needed. This mobility increases the flexibility of the workspace and makes the unit a valuable element in everyday clinical operations.