



Service unit 500X475X810h Handy-Soft, 4 drawers, Ral color handles

Available on request

SKU: DE9.31.4SO

Categories: Side Cabinets

Description

The **modular clinical side cabinet with 4 drawers**, measuring 500x475x810h, is designed to enhance organisation within clinical environments. Its compact size helps to optimise available space, offering order, practicality and reliable support for daily professional activities.

STRUCTURE AND FUNCTIONALITY

The cabinet features a clean and simple design, ideal for settings where hygiene and visual clarity are essential. The four drawers allow tools and materials to be arranged logically, ensuring they are always easy to access. This internal organisation improves workflow efficiency and reduces preparation times during clinical procedures.

CLINICAL MODULARITY

Thanks to its modular nature, the modular clinical side cabinet with 4 drawers integrates effortlessly into medical practices, aesthetic clinics and multidisciplinary environments. Its layout supports the operational needs of the clinician while contributing to a professional and reassuring atmosphere for patients.

MOBILITY AND OPERATIONAL SUPPORT

Equipped with wheels, the unit can be moved and positioned quickly according to procedural requirements. This mobility increases the flexibility of the workspace, making the cabinet a valuable support element in everyday clinical routines.