



Service unit 500X475X810h Handy-Soft, 5 drawers, Ral color handles

Available on request

SKU: DE9.31.5SO

Categories: Side Cabinets

Description

The **modular clinical side cabinet with 5 drawers**, measuring 500x475x810h, is designed to improve organisation within clinical environments. Its compact structure helps optimise space, providing order, practicality and reliable support throughout daily activities.

STRUCTURE AND FUNCTIONALITY

Featuring a clean and simple design, the cabinet suits settings where hygiene and clarity are essential. The five drawers allow a clear and efficient arrangement of instruments and frequently used materials, ensuring everything is within easy reach. This improves workflow and reduces preparation time during clinical procedures.

CLINICAL MODULARITY

Thanks to its modular configuration, the **modular clinical side cabinet with 5 drawers** integrates easily into a wide range of medical, aesthetic and multidisciplinary practices. Its internal organisation supports operational needs and contributes to a professional, tidy and reassuring clinical environment.

MOBILITY AND OPERATIONAL SUPPORT

Equipped with wheels, the unit can be moved quickly and positioned wherever needed, depending on procedural requirements. This mobility enhances workspace flexibility, making the cabinet an essential support element in daily clinical routines.