



## Service unit 500X475X810h Handy-Soft, 6 drawers, Ral color handles

Available on request

SKU: DE9.31.6SO

Categories: Side Cabinets

### Description

---

The **modular clinical side cabinet with 6 drawers**, measuring 500x475x810h, is designed to improve organisation within clinical environments. Its compact build allows for efficient use of available space, offering order, practicality and effective support throughout daily operations.

#### STRUCTURE AND FUNCTIONALITY

The cabinet features a clean, minimalist design, ideal for clinical settings where hygiene and visual clarity are essential. The six drawers provide well-divided storage areas, allowing instruments and frequently used materials to be arranged logically. The structure of the **modular clinical side cabinet with 6 drawers** makes workflows smoother, reduces preparation times and supports an efficient clinical routine.

#### CLINICAL MODULARITY

Thanks to its modular configuration, the unit integrates easily into various medical, aesthetic and multidisciplinary practices. Its internal organisation supports operational needs and contributes to a professional, tidy and reassuring environment for both clinicians and patients.

#### MOBILITY AND OPERATIONAL SUPPORT

Fitted with wheels, the cabinet can be repositioned quickly and placed wherever required. This mobility adds flexibility to the workspace, making the unit a valuable support element in day-to-day clinical procedures.