



Service unit 500X475X810h Handy-Soft, 6 drawers, Ral color handles

Available on request

SKU: DE9.31.6SO

Categories: Side Cabinets

Description

The **modular clinical side cabinet with 6 drawers**, measuring 500x475x810h, is designed to improve organisation within clinical environments. Its compact build allows for efficient use of available space, offering order, practicality and effective support throughout daily operations.

STRUCTURE AND FUNCTIONALITY

The cabinet features a clean, minimalist design, ideal for clinical settings where hygiene and visual clarity are essential. The six drawers provide well-divided storage areas, allowing instruments and frequently used materials to be arranged logically. The structure of the **modular clinical side cabinet with 6 drawers** makes workflows smoother, reduces preparation times and supports an efficient clinical routine.

CLINICAL MODULARITY

Thanks to its modular configuration, the unit integrates easily into various medical, aesthetic and multidisciplinary practices. Its internal organisation supports operational needs and contributes to a professional, tidy and reassuring environment for both clinicians and patients.

MOBILITY AND OPERATIONAL SUPPORT

Fitted with wheels, the cabinet can be repositioned quickly and placed wherever required. This mobility adds flexibility to the workspace, making the unit a valuable support element in day-to-day clinical procedures.