



Set of 2 Ivory SL4201C 20:1 green contra-angle with external spray

Kit consisting of two Ivory Silk Line 20:1 surgical contra-angles with external spray for implantology procedures in dental practice.

Available on request

SKU: OF420

Categories: Contra-Angles handpieces

Description

The Ivory Silk Line 20:1 surgical contra-angle kit with external spray includes two instruments intended for implantology procedures in dental practice.

The 20:1 reduction ratio delivers controlled speed and high torque during surgical use.

The precision-engineered transmission components are manufactured to very tight tolerances, ensuring efficient transfer of power from the micromotor while minimising vibration, noise and energy loss.

This results in smooth and stable performance throughout clinical procedures.

Compact head design improves access to difficult surgical areas, supporting better visibility and accuracy during treatment.

External spray system helps maintain effective cooling of the working area during use. T

The contra-angles are suitable for thermal disinfection and can be sterilised in an autoclave up to 135°C, in line with hygiene and sterilisation protocols in dental practice.

Technical features Ivory Silk Line 20:1 surgical contra-angle kit

Feature	Specification
Removable head	No
Fiber optics	No
Transmission ratio	20:1
Indications	Implantology



Feature	Specification
Locking mechanism	Push-button
Spray type	External
Compatibility	Intramatic ISO 3964
Rotation speed	2000rpm
Torque	? 55Ncm
Rotatable head	No
Head diameter	9.5mm
Head height	13.5mm
Bur diameter	2.334/2.350mm
Bur length	9–34mm
Weight	58g
Autoclavable	Yes (max 135°C)
Lubrication	Oil