



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/silver-seat-a6-daffodil-yellow-chrome-plated-base/>



SILVER seat, A6 daffodil yellow, chrome-plated base

SILVER seat, colour A6 daffodil yellow, chrome-plated base, height and backrest adjustment

364,74 €

Available on request

SKU: DE8.01.A6

Categories: Standard

Description

SILVER seat A6 daffodil yellow.

Comfortable and ergonomic, our Silver stools respect your body's anatomy and stand out for their excellent quality-price ratio. The height adjustment levers are easy to reach and let you adjust seat position without having to interrupt treatment.

Our range includes eight different colours for a perfect match with the finish and materials of your dental chair.

Soft, anatomic padding ensures maximum comfort even during long operations. Seamless, fire-repellent upholstery makes these stools easy to clean and disinfect.

All our stools are fitted with pneumatic height adjustment and backrest angle adjustment.

SILVER seat A6 daffodil yellow.