



## SILVER seat, A9 sapphire blue, chrome-plated base

SILVER seat, colour A9 sapphire blue, chrome-plated base, height and backrest adjustment

**Available**

**SKU:** DE8.01.A9

**Categories:** Standard

### Description

---

SILVER seat A9 sapphire blue.

Comfortable and ergonomic, our Silver stools respect your body's anatomy and stand out for their excellent quality-price ratio. The height adjustment levers are easy to reach and let you adjust seat position without having to interrupt treatment.

Our range includes eight different colours for a perfect match with the finish and materials of your dental chair.

Soft, anatomic padding ensures maximum comfort even during long operations. Seamless, fire-repellent upholstery makes these stools easy to clean and disinfect.

All our stools are fitted with pneumatic height adjustment and backrest angle adjustment.

SILVER seat A9 sapphire blue.