



## Single washbasin in stainless steel with surg. lever-operated mixer tap

Available on request

SKU: WASHCH-1

Categories: Sinks

### Description

---

The **stainless-steel single scrub sink** is designed to ensure **maximum hygiene and reliability** in operating theatres, clinics and sterilisation areas. Made entirely of **stainless steel**, it provides a **robust, corrosion-resistant and easy-to-clean structure**, ideal for sterile and high-demand medical environments.

The **single-station configuration** allows comfortable use during pre-surgical washing procedures, optimising ergonomics for medical personnel. Equipped with a **surgical lever mixer**, the sink enables **hands-free operation**, reducing the risk of cross-contamination and ensuring precise water control.

An **integrated soap dispenser** provides immediate access to cleaning agents, supporting proper surgical handwashing protocols.

The **deep bowl with rounded edges** ensures splash-free performance, while the **smooth, fully welded surfaces** simplify maintenance and prevent residue accumulation.

With its combination of **durability, functionality and ease of cleaning**, this stainless-steel single scrub sink offers a **professional-grade solution for hospitals, operating theatres and clinics** where hygiene and reliability are essential.