



## Stainless steel anaesthesia cart, 650x410x1150 mm

Shipping expected in approximately 21 days

SKU: TM5007

Categories: Carts

### Description

---

The **stainless-steel anaesthesia trolley** is designed to ensure **organisation, hygiene and safety** in operating theatres and clinical preparation areas. Built entirely from **stainless steel**, it provides long-term strength and resistance to corrosion, ideal for sterile environments requiring frequent cleaning.

Measuring **650x410x1150 mm**, the **stainless-steel anaesthesia trolley** features **six vertical sliding drawers**, a **lockable compartment with hinged door**, and an **open horizontal shelf** for frequently used materials. The **top surface is bordered on all four sides**, preventing instruments and containers from falling during use.

Mounted on **swivelling castors**, it moves smoothly and quietly even in confined spaces. Every detail has been developed to optimise workflow during anaesthetic and surgical procedures, while maintaining excellent **ergonomics and operational safety**.

The high-grade stainless-steel construction provides **antibacterial, easy-to-clean surfaces**, ensuring compliance with modern hospital hygiene standards. Perfect for **operating theatres, clinics and veterinary surgeries** requiring **durable, functional and professional-grade equipment**.