

## Stainless steel bag cart (no lid available)

Available on request

SKU: TM5042/1

Categories: Carts



### Description

---

The **stainless-steel sack trolley** is designed for **hygienic waste collection** in hospitals, clinics and laboratories.

Its **stainless-steel frame** ensures strength, easy cleaning and long service life.

**Smooth, seamless surfaces** prevent residue build-up, while **swivelling castors** allow quiet and effortless mobility.

Compact, ergonomic and durable, it is ideal for environments requiring **maximum hygiene and efficiency**.