



Stainless steel cart cabinet, dimensions 790x670xh140mm

Stainless steel cart cabinet, dimensions 790x670xh140mm, 2 internal shelves and antistatic casters

Shipping expected in approximately 21 days

SKU: TM2104

Categories: Carts

Description

The **stainless-steel container cabinet trolley** is designed for **efficient, hygienic and secure handling** of surgical containers and sterile materials. Made entirely of **stainless steel**, it offers exceptional strength, corrosion resistance and ease of cleaning—ideal for demanding medical environments.

Measuring **790x670xH1400 mm**, the trolley includes **two internal shelves** providing well-organised storage for surgical sets or sterilised instruments. The **smooth, fully welded surfaces** prevent dust and residue accumulation, ensuring quick and thorough sanitisation.

Fitted with **antistatic swivelling castors**, the **stainless-steel container cabinet trolley** moves **smoothly and quietly**, allowing precise positioning even in confined spaces. The closed-cabinet design protects the contents from dust and contamination, making it suitable for sterile departments, operating theatres and CSSD areas.

Combining ergonomic design and robust stainless-steel construction, this trolley provides a **reliable, hygienic and long-lasting solution** for clinics, laboratories and hospitals seeking professional-grade medical furniture.