



## Stainless steel operating table, dimensions 1500x600xh900mm (WxDxH)

Stainless steel veterinary operating table, dimensions  
1500x600xh900mm (WxDxH)

**Available on request**

**SKU:** TAO1500X

**Categories:** Surgical Operating Tables

### Description

---

The **stainless steel veterinary operating table 1500x600x900mm** provides a stable, hygienic and practical surface for daily clinical work. Made from 304 stainless steel, it offers excellent durability and allows quick disinfection, which is essential in sterile veterinary environments.

The table includes **swivelling wheels** that make movement easy and immediate. This mobility helps the clinician reposition the workstation quickly, improving workflow efficiency during examinations and procedures.

The 40mm worktop is built with 1.5mm stainless steel sheet and reinforced internally. This guarantees a stable surface even during delicate operative tasks. In addition, the recessed frame on all four sides allows the table to be positioned against a wall, even where coving is present. This helps maintain a well-organised and functional clinical space.

The **stainless steel veterinary operating table 1500x600x900mm** is suitable for clinics that need a durable, hygienic and easy-to-clean work surface. Its configuration supports both general examinations and more advanced veterinary procedures.