



Stainless steel operating table, dimensions 1800x600xh900mm (WxDxH)

Stainless steel veterinary operating table, dimensions 1800x600xh900mm (WxDxH)

Available on request

SKU: TAO1800X

Categories: Surgical Operating Tables

Description

Operating table stainless steel 1800x600xh900mm.

All of the veterinary examination tables are constructed from durable 304 grade stainless steel and are exceptionally easy to sanitise, making them ideal for sterile environments.

This selection of tables is designed to meet every veterinarians need and requirements for their practice.

40mm thick worktop, made with 1.5mm thick sheet metal and reinforced with stainless steel reinforcements.

Structure set back on all four sides to be able to go against the wall even in case of coving.