



Stainless steel operating table, dimensions 1800x600xh900mm (WxDxH)

Stainless steel veterinary operating table, dimensions
1800x600xh900mm (WxDxH)

Available on request

SKU: TAO1800X

Categories: Surgical Operating Tables

Description

The **stainless steel veterinary operating table 1800x600x900mm** provides a wide, hygienic and stable work surface for daily clinical procedures. Made from 304 stainless steel, it offers excellent durability and allows fast and effective disinfection, which is essential in sterile veterinary environments.

The 40mm worktop is built with 1.5mm stainless steel sheet and reinforced internally, ensuring stability even during delicate procedures. The recessed frame on all four sides allows the table to be placed directly against a wall, even where coving is present, helping maintain a tidy and functional clinical area.

The **stainless steel veterinary operating table 1800x600x900mm** is ideal for veterinary clinics that need a robust, hygienic and easy-to-clean surface. Its configuration supports both routine examinations and more advanced veterinary procedures.