

## Stainless steel operating table, dimensions 1800x600xh900mm (WxDxH)



Stainless steel veterinary operating table, dimensions 1800x600xh900mm (WxDxH)

**Available on request**

**SKU:** TAO1800X

**Categories:** Surgical Operating Tables

### Description

---

The **stainless steel veterinary operating table 1800x600x900mm** provides a wide, hygienic and stable work surface for daily clinical procedures. Made from 304 stainless steel, it offers excellent durability and allows fast and effective disinfection, which is essential in sterile veterinary environments.

The table features **swivelling wheels** that make repositioning easy and immediate. This mobility helps the clinician adapt the workstation quickly during examinations, preparations and operative tasks, improving workflow efficiency.

The 40mm worktop is built with 1.5mm stainless steel sheet and reinforced internally, ensuring stability even during delicate procedures. The recessed frame on all four sides allows the table to be placed directly against a wall, even where coving is present, helping maintain a tidy and functional clinical area.

The **stainless steel veterinary operating table 1800x600x900mm** is ideal for veterinary clinics that need a robust, hygienic and easy-to-clean surface. Its configuration supports both routine examinations and more advanced veterinary procedures.