



## Stainless steel operating table on casters, 1200x700xh900mm (WxDxH)

Stainless steel veterinary operating table on casters, dimensions 1200x700xh900mm (WxDxH)

**Shipping expected in approximately 21 days**

**SKU:** TAO1200XR

**Categories:** Surgical Operating Tables

### Description

---

The **stainless steel veterinary operating table 1200x700x900mm** is designed to provide a stable, hygienic and practical surface for daily clinical work. The 304 stainless steel structure ensures excellent durability and allows quick and thorough disinfection, which is essential in sterile veterinary environments.

The table is equipped with **swivelling wheels**, making movement easy and immediate. This mobility helps the clinician adjust the workstation quickly and improves workflow efficiency during examinations and procedures.

The worktop is 40 mm thick and made from 1.5 mm stainless steel sheet with reinforced support. This construction ensures stability during routine examinations and more demanding clinical tasks. In addition, the recessed frame on all four sides allows the table to be placed against a wall, even when coving is present. As a result, the clinical area becomes easier to organise and more efficient to manage.

The **stainless steel veterinary operating table 1200x700x900mm** is suitable for veterinary clinics and medical facilities that require a robust and easy-to-clean work surface. Its configuration supports both general examinations and more advanced operative procedures.