



Stainless steel operating table on casters, 1200x700xh900mm (WxDxH)

Stainless steel veterinary operating table on casters, dimensions 1200x700xh900mm (WxDxH)

Available on request

SKU: TAO1200XR

Categories: Surgical Operating Tables

Description

Operating table stainless steel 1200x700x900mm casters.

All of the veterinary examination tables are constructed from durable 304 grade stainless steel and are exceptionally easy to sanitise, making them ideal for sterile environments.

Mobile tables feature swivel castors for easy manoeuvrability.

This selection of tables is designed to meet every veterinarians need and requirements for their practice.

40mm thick worktop, made with 1.5mm thick sheet metal and reinforced with stainless steel reinforcements.

Structure set back on all four sides to be able to go against the wall even in case of coving.