



Stainless steel operating table on casters, 1500x700xh900mm (WxDxH)

Stainless steel veterinary operating table on castors, dimensions 1500x700xh900mm on castors (WxDxH)

Shipping expected in approximately 21 days

SKU: TAO1500XR

Categories: Surgical Operating Tables

Description

The **stainless steel veterinary operating table 1500x700x900mm** is designed to provide a stable, hygienic and practical surface for everyday clinical work. Made from 304 stainless steel, it offers excellent resistance and allows fast and complete disinfection, which is essential in sterile veterinary settings.

The table includes **swivelling wheels** that make movement simple and immediate. This mobility helps the clinician adapt the workstation quickly and improves workflow efficiency during procedures and examinations.

The 40mm worktop is built with 1.5mm stainless steel sheet and reinforced support. This construction ensures stability during routine visits and more complex clinical tasks. In addition, the recessed frame on all four sides allows the table to be placed directly against a wall, even where coving is present. As a result, the workspace becomes easier to organise and more functional.

The **stainless steel veterinary operating table 1500x700x900mm** is suitable for clinics that need a durable, easy-to-clean and professional operating surface. Its configuration supports both general examinations and advanced veterinary procedures.