



Stainless steel operating table on casters, 1500x700xh900mm (WxDxH)

Stainless steel veterinary operating table on castors, dimensions 1500x700xh900mm on castors (WxDxH)

Available on request

SKU: TAO1500XR

Categories: Surgical Operating Tables

Description

Operating table stainless steel 1500x700xh900mm casters.

All of the veterinary examination tables are constructed from durable 304 grade stainless steel and are exceptionally easy to sanitise, making them ideal for sterile environments.

This selection of tables is designed to meet every veterinarians need and requirements for their practice.

40mm thick worktop, made with 1.5mm thick sheet metal and reinforced with stainless steel reinforcements.

Structure set back on all four sides to be able to go against the wall even in case of coving.