



## Stainless steel washbasin cabinet,500x480X1020h,flush welded sink,basket

Available

SKU: WASHCST-5

Categories: Sinks

### Description

---

The **surgical stainless-steel wash unit** is designed for **hospitals, clinics and laboratories** requiring **maximum hygiene and durability**.

With dimensions **500x480x1020 H**, it combines a **compact structure** with high mechanical strength, ideal for confined or intensive-use environments.

The **welded flush-top bowl** ensures a **seamless, easy-to-clean surface**, eliminating crevices where dirt or bacteria could accumulate. The fully welded construction and rounded edges guarantee long-lasting stability and effortless maintenance.

A **foot-operated solenoid valve** allows **hands-free water control**, helping to maintain hygiene and efficiency during washing procedures.

An integrated **basket** provides organised storage for cleaning materials or accessories, improving workflow efficiency.

The **stainless-steel structure** resists corrosion and disinfectants, maintaining a high level of hygiene in **operating theatres and sterilisation rooms**.

Compact, functional and hygienic, this **surgical stainless-steel wash unit** represents a **professional solution** for healthcare environments that demand reliability, cleanliness and ergonomic design.