



Stainless steel washbasin cabinet, 500x480X900h, flush welded sink, basket

Available on request

SKU: WASHCH-ST

Categories: Sinks

Description

The **surgical stainless-steel wash unit 500x480x900 H** is designed for **medical, laboratory and clinical environments** that require **hygiene, functionality and durability**.

Made entirely of **stainless steel**, it ensures **robustness, corrosion resistance and easy maintenance**, even with frequent use.

The **welded flush-top bowl** provides a smooth, joint-free surface that prevents dirt accumulation and simplifies cleaning.

Rounded edges enhance operator safety and improve hygiene in everyday use.

A **foot-operated solenoid valve** allows **hands-free water control**, helping to maintain hygiene and efficiency during washing procedures.

The **integrated stainless-steel basket** offers convenient storage for cleaning materials and accessories, keeping the area organised.

Compact, hygienic and durable, this **surgical stainless-steel wash unit** is the ideal solution for **clinics, laboratories and medical rooms** that demand reliability, cleanliness and ergonomic design.