



## Stainless steel washbasin cabinet, 800X650X1020, with solenoid valve

Available on request

SKU: WASHCH-381

Categories: Sinks

### Description

---

The **surgical stainless-steel wash unit 800x650x1020 H** is designed for **hospitals, clinics and operating theatres** where hygiene and efficiency are essential.

Made entirely of **stainless steel**, it provides **robust construction, corrosion resistance and easy cleaning**.

The **welded flush-top bowl** creates a seamless surface that prevents dirt build-up and simplifies maintenance.

All corners are rounded and fully welded to enhance operator safety and ensure maximum hygiene.

A **foot-operated solenoid valve** allows **hands-free water control**, reducing the risk of cross-contamination.

An integrated **basket** provides convenient storage for cleaning accessories or consumables, keeping the workspace tidy.

Compact, reliable and easy to sanitise, this **surgical stainless-steel wash unit** delivers a **professional-grade solution** for healthcare facilities demanding ergonomic, hygienic and durable equipment.