



Stainless steel washbasin cabinet,800X650X1020, with solenoid valve

Available on request

SKU: WASHCH-381 **Categories:** Sinks

Description

The surgical stainless-steel wash unit 800×650×1020 H is designed for hospitals, clinics and operating theatres where hygiene and efficiency are essential.

Made entirely of stainless steel, it provides robust construction, corrosion resistance and easy cleaning.

The **welded flush-top bowl** creates a seamless surface that prevents dirt build-up and simplifies maintenance.

All corners are rounded and fully welded to enhance operator safety and ensure maximum hygiene.

A foot-operated solenoid valve allows hands-free water control, reducing the risk of cross-contamination.

An integrated **basket** provides convenient storage for cleaning accessories or consumables, keeping the workspace tidy.

Compact, reliable and easy to sanitise, this **surgical stainless-steel wash unit** delivers a **professional-grade solution** for healthcare facilities demanding ergonomic, hygienic and durable equipment.