


## Stand for 2 Pyrex glasses for CLEA3L, CLEA5L, size 262x162



Holder for 2 Pyrex glasses for CLEA3L CLEA5L, 8 mm thick, size 262x162

**Available**

**SKU:** CLA04

**Categories:** Beakers

### Description

---

Stand for 2 glasses.

Stand for 2 Pyrex glasses for Ultrasonic Cleaners Digiclean CLEA3L (3lt.) and CLEA5L (5lt.).

Dimensions 262 x 162 mm.