


Sterilization labeller with 500 label roll

Sterilization labeller, 2 print lines, steam indicator, with 500 label roll
26x16mm

Available

SKU: DEA4.13

Categories: Packaging Equipment



Description

Printex Medical labeler has two print lines and is the ideal appliance for tracking sterilized products.

The labels apply after packaging as they are equipped with an indicator that identify the sterilization.

Featuring a double adhesive support the label can simply be attached to the sterile wrap and applied to patient documentation.

Maximum security thanks to the ability to enter the sterilization date, the expiration date and the lot number.

Printex Medical labeler for sterilization, with 2 printing lines, Roll with 500 labels included.