

Sterilization labeller with 500 label roll

Sterilization labeller, 2 print lines, steam indicator, with 500 label roll 26x16mm

Available from 18/05/2026

SKU: DEA4.13

Categories: Labeling Machines



01.09.14	890	01.09.14	890
30.09.14	890	30.09.14	890

Description

Sterilization Labeler – Printex Medical, 2-Line Print Head with 500-Label Roll Included

The **Printex Medical sterilization labeler** with a 2-line print head is the ideal solution for traceability of sterilized instruments in medical and dental settings. Designed for clinical environments that demand precision and compliance, this tool ensures reliable labeling throughout the sterilization process.

Each label features a **chemical indicator** that confirms successful sterilization by changing color. This provides an immediate visual check, improving safety and workflow efficiency.

Thanks to the **dual-adhesive backing**, labels can be easily removed from sterile packaging and reapplied directly to the patient's medical records. This supports proper documentation and long-term traceability—critical for regulatory compliance and risk management.

The two-line print layout allows users to add key data such as the **sterilization date, expiration date, and batch number**—ensuring every step is documented with accuracy and clarity.

Built for heavy daily use, the labeler combines compact design with durability. It comes equipped with a **roll of 500 high-adhesion sterilization labels**, making it ready to use straight out of the box.

Whether for a dental practice, outpatient clinic, or surgical center, the Printex Medical labeler is a must-have tool for safe, organized, and compliant sterilization tracking.