




Printed on 30/05/2025

<https://shop.tecnomeditalia.com/en/product/sterilization-roll-10-cm-high-200-m-long-pack-4-pieces/>



Sterilization roll 10 cm high, 200 m long. (Pack 4 pieces)

Sterilization roll 10 cm high and 200 m long. Pack of 4.

105,02 €

Available on request

SKU: DEA4.11

Categories: Consumable

Description

Sterilization roll 10 cm high, 200 m long. (Pack 4 pieces)

Tecnomed Italia offers a range of sterilization rolls that ensure a high level of protection against bacteria.

High quality materials: the rolls are composed of Medical grade paper from 60 to 70 g/square meter and a reinforced laminated film that guarantee an excellent opening phase without lacerations of the plastic material or possible release of paper fibers.

Triple welding line: the constant sealing during sterilization is ensured by the triple welding line.

Appropriate packaging: the rolls are individually packaged in protective plastic film and placed into carton boxes. This ensures dust protection even after the carton box is opened

Class 1 Process indicators (UNI EN ISO 11140-1): Indicators change color when sterilization has occurred

Compliant with the rules: all materials comply with international standards.

Features Sterilization roll 10 cm.

- Transparent film rolled out of blue color that allows easy identification of the content.
- Reinforced film to avoid lacerations of plastic during opening
- 60/70 gsq Medical grade paper of superior quality.
- The triple welding line is characterized by resistance and spacing both ensuring constant tightness during sterilization process as well as excellent peelable during opening without release of paper fibers or lacerations.
- Process indicators for steam and EO
- Non-toxic water inks certified ISO 11140
- All materials are in compliance with international standards
- Indication of opening direction

- Validated production process

The rolls are individually packaged in protective plastic film and placed into carton boxes. This ensures dust protection even after the carton box is opened.

Suitable for all sealing machines.