



## STING ST seat, 4009 cyclamen, aluminium base

STING ST seat, colour 4009 cyclamen, aluminium base, height adjustment and backrest rotation

**Available on request**

**SKU:** DE8.05T3492

**Categories:** Over Line

### Description

---

STING ST seat 4009 cyclamen.

An appropriate stool is an important investment for the comfort and health of the dental professional. Tecnomed has a vast experience in the design and production of ergonomic dental stools. Tecnomed stools help you maintain a correct posture with natural curvature of the spine, uniform distribution of body weight and no painful compressions.

StingST stools are designed to offer superior comfort for dentists and patients at all times.

StingST stools are upholstered in top quality imitation leather and are free from stitching for maximum hygiene and easy cleaning.

The upholstery material is fireproof according to ISO 6875 and UNI 4818.74. This top quality upholstery helps ensure a soft, anatomical seating position with uniform weight distribution and no painful compressions.

The StingST stool is specially designed for dentists who work in a seated position for many hours at a time. Its special saddle shape seat promotes correct posture. The “V” leg position also relieves strain on the spine and improves circulation. Seat angle is adjustable to suit requirements.

STING ST seat 4009 cyclamen.