



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/stool-beauty-steel-basement-seatregulation-tecnomed/>



Stool, Beauty, steel basement, seatregulation, Tecnomed

Beauty Stool Tecnomed, steel basement, with seat regulation, ergonomic structure, available in 30 colors.

646,27 €

Available

SKU: DE8.07X

Categories: Handy Line

Description

Beauty stool.

An appropriate stool is an important investment for the comfort and health of the dental professional.

Beauty is a professional dentist chair designed for complete comfort even after many hours of work.

Its ergonomic structure helps you to maintain correct posture while distributing weight in a uniform manner.

Their antibacterial and water-repellent imitation leather upholstery makes these chairs extremely easy to disinfect.

The range offers a choice of thirty different colours for a perfect match with the finish and materials of your dental chair.

Top quality imitation leather upholstery ensures a soft, anatomical seating position with uniform weight distribution and no painful compressions.

Beauty stool chair feature adjustable pneumatic height adjustment and backrest depth adjustment.