


Support vesa 75/100 dental white



VESA fixing plate for connecting monitors with 70x70mm or 100x100mm interface to monitor support arms using M5 screw fixing.

Available

SKU: SV029010

Categories: VESA Mounts

Description

SV029010 – VESA mounting plate for monitor support arm.

Main features

- VESA interface 70x70mm and 100x100mm
- Clamping with 2 x M5 screws

Technical specifications

- Compatible VESA standards: 70x70mm / 100x100mm
- Material: metal
- Colour: dental white
- Certification: 2001/95/EC
- Warranty: 36 months

Applications

- Dental / podiatry units
- Medical workstations with monitor support arm

Configuration note

The VESA mounting plate for monitor support arm can be combined with **Tecnomed Italia monitor support arms**.