




Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/support-vesa-75-100-dental-white/>

Support vesa 75/100 dental white



VESA 75/100 mouth support for monitor arm mounting, colour dental white, dimensions (LxH) 111x115 mm

35,34 €

Available

SKU: SV029010

Categories: Professional Equipment

Description

Support VESA 75/100.

Though they are indispensable items of equipment in any modern dental studio, monitors and computers must not compromise the high levels of hygiene needed. Medical PCs and monitors must therefore conform to strict hygiene standards as well as providing a high level of functionality in terms of mobility and usability.

Stylish, attractive looks are also a plus. Our arms support all types of monitor. Installation is extremely easy and does not require specialist assistance, just a simple Allen key! They also have a convenient handle so that you can move the monitor or all in one PC around easily.

Fixing holes are provided for Vesa 75 and Vesa 100 mountings. This new range of modular monitor arms from Tecnomed Italia offers amazing versatility and lets you specify components of various sizes. You can therefore configure your monitor arm exactly as you want it. Cables can be routed through the arm as an option. Various cables are available, including VGA, USB and power

Technical features Support VESA 75/100.

Color dental white

L 115 mm

H 115 mm

Warranty 24 months