


## Support vesa 75/100 dental white



VESA fixing plate for connecting monitors with 70×70mm or 100×100mm interface to monitor support arms using M5 screw fixing.

**Available**

**SKU:** SV029010

**Categories:** VESA Mounts

### Description

---

**SV029010 – VESA mounting plate for monitor support arm.**

#### Main features

- VESA interface 70×70mm and 100×100mm
- Clamping with 2 × M5 screws

#### Technical specifications

- Compatible VESA standards: 70×70mm / 100×100mm
- Material: metal
- Colour: dental white
- Certification: 2001/95/EC
- Warranty: 36 months

#### Applications

- Dental / podiatry units
- Medical workstations with monitor support arm

#### Configuration note

The VESA mounting plate for monitor support arm can be combined with **Tecnomed Italia monitor support arms**.