

Support vesa 75/100 with handle dental white



VESA 75 /100 mounting plate with integrated handling handle for connection of the monitor to support arms.

Available

SKU: SV019010

Categories: VESA Mounts

Description

SV019010: VESA mounting plate with handling handle for monitor support arm.

Main features

- VESA mounting pattern 75x75mm and 100x100mm
- Integrated handle for monitor positioning and movement
- Fastening with 2 x M5 screws

Technical specifications

- Compatible with VESA standards 75x75mm / 100x100mm
- Handle dimensions (LxH): 290mm x 170mm
- Material: metal
- Colour: dental white
- Certification: 2001/95/EC



- Warranty: 36 months

Applications

- Dental / podiatry units
- Medical workstations with articulated arm

Configuration note

The VESA mounting plate with handle can be combined with Tecnomed monitor support arms.