

Support vesa 75/100 with handle dental white



VESA support 75/100 with handle for monitor arm mounting, color dental white

Available

SKU: SV019010

Categories: VESA Mounts

Description

VESA support 75/100 with handle

Though they are indispensable items of equipment in any modern dental studio, monitors and computers must not compromise the high levels of hygiene needed. Medical PCs and monitors must therefore conform to strict hygiene standards as well as providing a high level of functionality in terms of mobility and usability.

Stylish, attractive looks are also a plus. Our arms support all types of monitor. Installation is extremely easy and does not require specialist assistance, just a simple Allen key! They also have a convenient handle so that you can move the monitor or all in one PC around easily.

Fixing holes are provided for Vesa 75 and Vesa 100 mountings. This new range of modular monitor arms from Tecnomed Italia offers amazing versatility and lets you specify components of various sizes. You can therefore configure your monitor arm exactly as you want it. Cables can be routed through the arm as an option. Various cables are available, including VGA, USB and power

Technical features VESA support 75/100 with handle

Color	dental white
Handle	•
A	115 mm
B	115 mm
L	245 mm
Warranty	24 months

