



Printed on 01/06/2025

<https://shop.tecnomeditalia.com/en/product/supporto-a-parete-per-braccio-okkio-colore-bianco/>



Wall support for OKKIO arm, white colour

<p class="p1">Support for wall mounting of OKKIO arm, white colour.</p>

205,31 €

Available

SKU: FK934P

Categories: Arm for cameras

Description

Support for wall mounting of OKKIO arm, white colour

For wall mounting, the best solution for stable video recording.

- maximum distance from Dental Chair headrest- 1,5 mt;
- white colour