




Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/supporto-a-poggiatesta-per-braccio-okkio/>



## Headrest support for OKKIO arm

Headrest support for OKKIO arm, for dental chair mounting, 180° degree rotatable support, minimum space required.

**247,28 €**

**Available**

**SKU:** FK939S

**Categories:** Arm for cameras

### Description

---

Headrest support for OKKIO arm, minimum space required, 180° degree rotatable support.

For dental chair mounting, a solution that ensures stable video recording.