



Supragingival powder tank for AP-A,AP-B DTE-Woodpecker, 65 UM

Supragingival powder tank for AP-A air flow air polisher, and AP-B ultrasonic scale and air polisher DTE Woodpecker, 65 UM.

Estimated shipping in 30 days

SKU: PL40APB65

Categories: Accessories

Description

Supragingival powder tank for AP-A and AP-B DTE.

Supragingival powder tank for AP-A air flow air polisher, and AP-B ultrasonic scale and air polisher DTE Woodpecker, 65 UM.