



Printed on 01/06/2025

<https://shop.tecnomeditalia.com/en/product/surface-disinfectant-metasys-mk-start-2/>

Surface disinfectant Metasys MK, start 2

Surface disinfectant Metasys MK, start 2, 1 dispenser with 70 wipes plus 6x70 wipes and a 5-litre canister SK

180,08 €

Available on request

SKU: 60030131

Categories: Surfaces



Description

Surface disinfectant Metasys MK start 2.

Alcohol-free disinfectant wipes, ready to use, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- **Broad-spectrum action:**
 - **HBV, HIV, HCV:** disinfects in just 1 minute.
 - **Bactericidal and fungicidal:** effective in 1 minute.

Certification: CE 0044.

Product Code: 60030130

Packaging Surface disinfectant Metasys MK start 2: 1 dispenser with 70 wipes plus 6x70 wipes and a 5-litre canister SK.