



Surface disinfectant Metasys MK, start 2

Surface disinfectant Metasys MK, start 2, 1 dispenser with 70 wipes plus 6x70 wipes and a 5-litre canister SK

Available on request

SKU: 60030131

Categories: Surface Disinfection

Description

Surface disinfectant Metasys MK start 2.

Alcohol-free disinfectant wipes, ready to use, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

Key Features:

- Long-lasting effect: prevents the settlement of new airborne germs.
- Advanced formula: new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- Broad-spectrum action:
 - HBV, HIV, HCV: disinfects in just 1 minute.
 - o Bactericidal and fungicidal: effective in 1 minute.

Certification: CE 0044.

Product Code: 60030130

Packaging Surface disinfectant Metasys MK start 2: 1 dispenser with 70 wipes plus 6×70 wipes and a 5-litre canister SK.