



Printed on 03/06/2025

<https://shop.tecnomeditalia.com/en/product/surface-disinfectant-metasys-mt-kit-2/>

## Surface disinfectant Metasys MT, kit 2

Metasys MT surface disinfectant, kit 2, 1 dispenser with 70 wipes plus 5x70 wipes and a 5-litre SD canister

**193,98 €**

**Available on request**

**SKU:** 60030094

**Categories:** Surfaces



### Description

---

Metasys MT surface disinfectant kit 2

**Alcohol-free disinfectant wipes**, ready to use, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

#### Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- **Broad-spectrum action:**
  - **HBV, HIV, HCV:** disinfects in just 1 minute.
  - **Bactericidal and fungicidal:** effective in 1 minute.

**Certification:** CE 0044.

**Product Code:** 60030130

**Packaging Metasys MT surface disinfectant kit 2:** Extra large wipes (290 x 270 mm), for use with GREEN&CLEAN SD, for wipe disinfection of large surfaces of medical devices.