

Surface disinfectant Metasys MT, kit 2

Metasys MT surface disinfectant, kit 2, 1 dispenser with 70 wipes plus 5x70 wipes and a 5-litre SD canister

Estimated shipping in 15 days

SKU: 60030094

Categories: Surface Disinfection



Description

Metasys MT surface disinfectant kit 2

Alcohol-free disinfectant wipes, ready to use, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- **Broad-spectrum action:**
 - **HBV, HIV, HCV:** disinfects in just 1 minute.
 - **Bactericidal and fungicidal:** effective in 1 minute.

Certification: CE 0044.

Product Code: 60030130

Packaging Metasys MT surface disinfectant kit 2: Extra large wipes (290 x 270 mm), for use with GREEN&CLEAN SD, for wipe disinfection of large surfaces of medical devices.