



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/surface-disinfectant-metasys-mt-refill-2/>

Surface disinfectant Metasys MT, refill 2

Surface disinfectant Metasys MT, refill 2, 1 dispenser with 70 wipes and 5 x 70 wipes

117,01 €

Available from 12/06/2025

SKU: 60030092

Categories: Surfaces



Description

Surface disinfectant Metasys MT refill 2

Alcohol-free disinfectant wipes, ready to use, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- **Broad-spectrum action:**
 - **HBV, HIV, HCV:** disinfects in just 1 minute.
 - **Bactericidal and fungicidal:** effective in 1 minute.

Certification: CE 0044.

Packaging Surface disinfectant Metasys MT refill 2: dispenser with 70 Extra large wipes (290 x 270 mm) + 750 ml SK bottle (MK START).