



## Surface disinfectant Metasys MT, refill 2

Surface disinfectant Metasys MT, refill 2, 1 dispenser with 70 wipes and 5 x 70 wipes

## **Available**

SKU: 60030092

Categories: Surface Disinfection

## **Description**

Surface disinfectant Metasys MT refill 2

**Alcohol-free disinfectant wipes**, ready to use, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

## **Key Features:**

- Long-lasting effect: prevents the settlement of new airborne germs.
- Advanced formula: new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- Broad-spectrum action:
  - o HBV, HIV, HCV: disinfects in just 1 minute.
  - o Bactericidal and fungicidal: effective in 1 minute.

Certification: CE 0044.

Packaging Surface disinfectant Metasys MT refill 2: dispenser with 70 Extra large wipes (290 x 270 mm) + 750 ml SK bottle (MK START).