

## Surface disinfectant Metasys MT, refill 2



Surface disinfectant Metasys MT, refill 2, 1 dispenser with 70 wipes and 5 x 70 wipes

**Available on request**

**SKU:** 60030092

**Categories:** Surfaces

### Description

---

Surface disinfectant Metasys MT refill 2

**Alcohol-free disinfectant wipes**, ready to use, designed for the disinfection of alcohol-sensitive surfaces on medical devices, such as dental units, handles, and leather/synthetic leather upholstery.

#### Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of quaternary ammonium compounds, free from aldehydes, phenols, and phosphates.
- **Broad-spectrum action:**
  - **HBV, HIV, HCV:** disinfects in just 1 minute.
  - **Bactericidal and fungicidal:** effective in 1 minute.

**Certification:** CE 0044.

**Packaging Surface disinfectant Metasys MT refill 2:** dispenser with 70 Extra large wipes (290 x 270 mm) + 750 ml SK bottle (MK START).