



Surface disinfectant Metasys SD, refill 1

Surface disinfectant Metasys SD, refill 1, 6 pouches x 750 ml

Available on request

SKU: 60030002

Categories: Surface Disinfection

Description

Surface disinfectant Metasys SD refill 1.

Alcohol-based nebulizing disinfectant, ideal for the disinfection of alcohol-resistant surfaces on medical devices, such as dental units, tables, cabinets, and receptacles.

Key Features:

- Long-lasting effect: prevents the settlement of new airborne germs.
- Advanced formula: new combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.
- Broad-spectrum action:
 - o HBV, HIV, HCV: disinfects in just 30 seconds.
 - o Bactericidal and fungicidal: effective in 1 minute.
 - o Tuberculocidal: action within 5 minutes.

Certification: CE 0044.

Packaging Surface disinfectant Metasys SD refill 1: 6 pouches x 750 ml.