



Printed on 02/06/2025

<https://shop.tecnomeditalia.com/en/product/surface-disinfectant-metasys-sd-refill-1/>

Surface disinfectant Metasys SD, refill 1

Surface disinfectant Metasys SD, refill 1, 6 pouches x 750 ml

88,70 €

Available on request

SKU: 60030002

Categories: Surfaces



Description

Surface disinfectant Metasys SD refill 1.

Alcohol-based nebulizing disinfectant, ideal for the disinfection of alcohol-resistant surfaces on medical devices, such as dental units, tables, cabinets, and receptacles.

Key Features:

- **Long-lasting effect:** prevents the settlement of new airborne germs.
- **Advanced formula:** new combination of alcohols, quaternary ammonium compounds, and glyoxal, free from formaldehyde.
- **Broad-spectrum action:**
 - **HBV, HIV, HCV:** disinfects in just 30 seconds.
 - **Bactericidal and fungicidal:** effective in 1 minute.
 - **Tuberculocidal:** action within 5 minutes.

Certification: CE 0044.

Packaging Surface disinfectant Metasys SD refill 1: 6 pouches x 750 ml.